


FineCross® MG


Coronary Micro-Guide Catheter

Finecross® MG microcatheter represents the integration of excellent crossability and optimal guidewire support to treat very challenging cases, involving tortuous vessels, CTO, bifurcations, highly stenosed and calcified after lesions*.



PRODUCT SPECIFICATIONS

Shaft coating:	Hydrophilic coating
Braiding:	Stainless steel - SUS
Inner layer:	Polytetrafluoroethylene - PTFE
Radiopaque marker:	0.7 mm gold marker located at 0.7 mm from the tip
Usable length:	130 or 150 cm




* Bench mark test result from Ashitaka factory




FineCross® MG

Coronary Micro-Guide Catheter



Conventional microcatheter



Finecross MG
Tapers from 2.6 Fr to 1.8 Fr.
over entire catheter length

ORDERING INFORMATION

	Distal outer diameter	Proximal outer diameter	Distal inner diameter	Proximal inner diameter	Catheter length	Guidewire compatibility
NC-F863A					130 cm	
NC-F865A	1.8 Fr / 0.60 mm	2.6 Fr / 0.87 mm	0.018" / 0.45 mm	0.021" / 0.55 mm	150 cm	Maximum diameter 0.014" / 0.36 mm